



Item # _____

Job _____



Starsys™ System Secure IV Solution Cart

STARSYS™ SYSTEM SECURE IV CART

The Starsys Secure IV Cart is a system of polymer components that offers security in a flexible storage solution.

Comprehensive System: Two different sized carts are constructed from a combination of polymer components.

Increases Productivity: Removable totes allow for the creation of an inventory replacement system with prestocked totes. Polymer construction prevents damage to products which results in additional work and cost.

Easily Accessible: Full extension shelves allow for solution bags to be stored standing up or laying down while providing complete access to all of the product.

JCAHO Compliant: The locking feature is in compliance with JCAHO's medication storage standards.

Visible Security: A locking clear door provides for visual inventory while contents are secure.

Enclosed Protection: Solid enclosure panels and a door protects contents from splashes or other potential contamination.

Corrosion Resistant: Constructed of durable corrosion-proof polymer panels that are easy to clean won't chip, rust, or dent.

Easily Maneuverable: The 58" height lighter polymer design and 5" (127mm) casters, allow for easier, safer transport.



SXSIVSTOR



InterMetro Industries Corporation
North Washington Street
Wilkes-Barre, PA 18705
www.metro.com



Job _____



Starsys™ System Secure IV Solution Cart

STARSYS™ SYSTEM SECURE IV CART

Dimensions

Cat. No.	Description	Height		Width		Length		Total Totes Per Cart
		(in.)	(mm)	(in.)	(mm)	(in.)	(mm)	
SXSIVSTOR	Single-Wide	58	1473	27 7/8	708	22 11/16	576	3
SXDIVSTOR	Double-Wide	58	1473	27 7/8	708	42 1/2	1080	6

- Units feature four 5" (127 mm) casters - the front two casters have color-matched toe-brakes.
- The dimensions of the tote boxes (MTB93080W) that are included are 8 x 17 1/2 x 22 1/2 inches (203 x 445 x 572 mm).

Accessories

Dividers

Cat. No.	Description	For Tote Box
MDS93080NAT	Short Dividers	MTB93080W
MDL93080NAT	Long Dividers	MTB93080W



Additional/Exchange Totes

Cat. No.	Description	Outside Dimensions (includes lip) Length x Width		Inside Dimensions Length x Width		Height O.D. (in.) (mm)	Capacity		Approx. Tote Box Wt. (lbs.)
		(in.)	(mm)	(in.)	(mm)		(cu. ft.)	(cu.m.)	
MTB93080W	8" White Tote Box	22 1/2 x 17 1/2	572x455	20 1/4 x 15 1/8	514x384	8 203	1.32 0.037	5.25	

Additional Shelves

Cat. No.	Description	Overall Dimensions Height/Depth/Width		Max. Wt. Capacity Uniformly Distributed	
		(in.)	(mm)	(in.)	(mm)
SXSWAL	Full Extension Single-Wide Shelf-LH	1 1/2 x 21 9/16 x 16 1/2	38 x 548 x 419	60	27.2
SXSWAR	Full Extension Single-Wide Shelf-RH	2 1/2 x 21 9/16 x 16 1/2	38 x 548 x 420	60	27.2
SXDWAL	Full Extension Double-Wide Shelf-LH	1 1/2 x 21 9/16 x 36 1/2	38 x 548 x 924	125	56.6
SXDWAR	Full Extension Double-Wide Shelf-LR	2 1/2 x 21 9/16 x 36 1/2	38 x 548 x 924	125	56.6

Snap-On Card Holder

Cat. No.	Description	(in.)	(mm)
OP2501CLR	Snap-On Card Holder	4 1/2 x 8	114 x 199

Note:

Other configurations can be built-to-order.
Contact your Metro rep for assistance.

Manufactured by:



InterMetro Industries Corporation
North Washington Street, Wilkes-Barre, PA 18705
Phone: 570-825-2741 • Fax: 570-825-2852
For Product Information Call: 1-800-433-2232
Visit Our Web Site: www.metro.com

Information and specifications are subject to change without notice. Please confirm at time of order.
Copyright © 2006 InterMetro Industries Corp.

EHC6-06

