



# STARSYS™



Overhead cabinet shelves are removable, easy-to-clean, and are available in both solid polymer and epoxy-coated, open-wire shelf designs.

Starsys mobile workcenter units can be specified with or without overhead cabinets attached. Available overhead heights include 24", 27", 30", and 36".

Overhead cabinet accessories include light fixture, cassette bins for small item storage and organization, and poly shelf dividers.

Overhead cabinets can be specified with a variety of available doors, clear or solid style, locking or non-locking.

Color-coded drawer pulls are available in an array of choices to fit any system or decor.

Standard countertop choices include laminate, stainless steel, phenolic and solid surface.

Starsys Doors, available in multiple heights, can be specified with or without window, as well as locking and non-locking.

Microban® antimicrobial protection inhibits the growth of stain and odor-causing bacteria on Starsys, keeping the surface areas "cleaner between cleanings."

Additional caster options are available through Starsys configurator.com or by contacting your Metro representative.

Starsys drawer options include choice of polymer drawers or heavy-duty drawers featuring stainless steel drawer interior. Heavy-duty drawers are available in both single- and double-wide modules.

5" (127mm) polymer total-lock casters provide a stable work station when locked and smooth rolling cart when not.



Starsys Mobile WorkCenters are built for change. They allow you to easily adapt to a changing work environment and even inspire change in your work environment.

**Total Unit Height:** All units with overhead storage have a total height of 82" (2083mm) (with 5" casters). Deduct 2" from unit height when specifying 3" Total-Lock casters or Leveling casters.

**Available Widths:** Starsys Mobile WorkCenters are available in four widths; single: 24" (610mm), double: 44" (1118mm), triple: 63" (1600mm) & single-kneewell-single: 73" (1854mm).

**Working heights with Overheads:** 33", 36", 39", 42" (45" height is not available with overhead storage)

**Kneewell Options:** Starsys Mobile WorkCenters configurations have several kneewell options including a keyboard tray, pencil drawer and support bracket. The triple unit has a right oriented kneewell option with a nominal opening of 42" (1067mm) and the single-kneewell-single configuration has a nominal opening of 30" (762mm).

**Working Heights:** Available working heights without overhead storage: 33", 36", 39", 42", 45" (with 5" casters). Deduct 2" (51mm) from unit height when specifying 3" Total-Lock casters or Leveling casters.

<p><b>3" Total-Lock Caster:</b> Utilize 3" total-lock plate casters for limited mobile applications or when additional storage area is required.</p>  <p>B3P-TL</p>	<p><b>5" Total-Lock, All-Polymer Caster:</b> Recommended for corrosive environments. All-polymer total-lock plate caster is an economical alternative to stainless casters.</p>  <p>B5PC-TL</p>	<p><b>Stabilizer/Leveling Caster:</b> Recommended for applications where benchtop equipment demands unit stability and/or a level worksurface. Engage caster foot for stability, release for mobility.</p>  <p>B3N-AS</p>
--	--	--

Starsys Carts are built to exact customer specifications, therefore, they are non-cancelable, non-returnable.

Preconfigured Mobile WorkCenters  
Mobile WorkCenters with Laminate Tops



Cat. No. **SXRSENT1**



Cat. No. **SXRDEMENT2**



Cat. No. **SXR40MW10**



Cat. No. **SXRT38MW14**

Alternate top materials and finishes available.



Cat. No. **SXRK32MW5**



Cat. No. **SXRK32MW6**



Cat. No. **SXRK41MW7**

To configure a Starsys unit to your exact specifications, please visit [www.metroconfigurator.com](http://www.metroconfigurator.com)



**SXRSENT1**

**SXRSENT1 includes:**

Single-Wide Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 Laminate Countertop with 3" (76mm) Backsplash (Wilsonart Carbon EV)  
Mobile Bay (Keylock/2 keys)  
 (2) 3" (76mm) Polymer Drawers  
 (2) 6" (152mm) Polymer Drawers  
 (1) 9" (229mm) Polymer Drawer  
Overhead Cabinet (27" (686mm) Interior Height)  
 (2) Polymer Shelves  
 Left-Hinged 225° Clear Locking Door (Keylock/2 keys)  
 Slate Blue Drawer & Door Pulls  
 Overhead Uprights & Crossbars  
Countertop Height: 39.5"H (1003mm)  
Countertop to Overhead Clearance: 16" (406mm)  
Overall Dimensions: 81.5"H (2070mm) x 23.6"W (600mm) x 28.6"D (726mm)



**SXRDEMENT2**

**SXRDEMENT2 includes:**

Double-Wide Locking Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, (2) Swivel)  
 Laminate Countertop with 3" (76mm) Backsplash (Wilsonart Carbon EV)  
Left Mobile Bay (Keylock/2 keys)  
 (2) 3" (76mm) Polymer Drawers  
 (2) 6" (152mm) Polymer Drawers  
 (1) 9" (229mm) Polymer Drawer  
Right Mobile Bay  
 (3) Polymer Shelves  
 Right Hinged 225° Solid Locking Door (Keylock/2 keys)  
Double-Wide Overhead Cabinet (27" (686mm) Interior Height)  
 (2) Double Wide Poly Shelves  
 Center-Closing 225° Clear Locking Door Set (Keylock/2 keys)  
 Slate Blue Drawer & Door Pulls  
 Overhead Uprights & Cross Bars  
Countertop Height: 39.5"H (1003mm)  
Countertop to Overhead Clearance: 16" (406mm)  
Overall Dimensions: 81.5"H (2070mm) x 43.4"W (1103mm) x 28.6"D (727mm)



**SXR40MW10**

**SXR40MW10** includes:

- Double-Wide Locking Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)
- Laminate Countertop with 3" (76mm) Backsplash (Wilsonart Carbon EV)
- Left Mobile Bay (Keylock/2 keys)
- (3) 6" (152mm) Polymer Drawers
- (1) 9" (229mm) Polymer Drawer
- Right Mobile Bay (Keylock/2 keys)
- (3) 6" (152mm) Polymer Drawers
- (1) 9" (229mm) Polymer Drawer
- Double-Wide Overhead Cabinet (27" (686mm) Interior Height)
- (2) Double-wide Polymer Shelves
- Center-Closing 225° Clear Locking Door Set (Keylock/2 keys)
- Slate Blue Drawer & Door Pulls
- Overhead Uprights & Crossbars
- Countertop Height: 39.5"H (1003mm)
- Countertop to Overhead Clearance: 16" (406mm)
- Overall Dimensions: 81.5"H (2070mm) x 43.4"W (1102mm) x 28.6"D (726mm)



**SXRT38MW14**

**SXRT38MW14** includes:

- Triple-Wide Non-Locking Mobile WorkCenter with 3" (76mm) Total Lock Casters
- Laminate Countertop with 3" Backsplash (Wilsonart Carbon EV)
- Left Mobile Bay
- (1) 3" (76mm) Polymer Drawer
- (3) 6" (152mm) Polymer Drawers
- (1) 9" (229mm) Polymer Drawer
- Right Mobile Bay
- (2) Double-wide Full –Extension Chrome Wire Shelves
- Center-closing 225° Clear Non-Locking Door Set
- Triple-Wide Overhead Cabinet (27" (686mm) Interior Height)
- Left Bay
- (2) Single-Wide dark taupe epoxy-coated Wire Shelves
- Left-Hinged 225° Clear Non-locking Door
- Right Overhead Bay
- (2) Double-Wide dark taupe epoxy-coated Wire Shelves
- Center-Closing 225° Clear Non-locking Door Set
- Slate Blue Drawer & Door Pulls
- Overhead Uprights & Crossbars
- Countertop Height: 37.5"H (953mm)
- Countertop to Overhead Clearance: 16" (406mm)
- Overall Dimensions: 81.5"H (2070mm) x 63.2"W (1606mm) x 28.6"D (726mm)



**SXRK32MW5**

**SXRK32MW5** includes:

73" (1855mm) Mobile WorkCenter (Non-Locking) with 27" (686mm) Knee-hole & 3" (76mm) Total Lock Casters  
 Laminate Countertop with 3" (76mm) Backsplash (Wilsonart Carbon EV)  
 Pencil Drawer & Knee-hole Connecting Brackets

Left Mobile Bay

- (1) 3" (76mm) Polymer Drawer
- (2) 6" (152mm) Polymer Drawers
- (1) 9" (229mm) Polymer Drawer

Right Mobile Bay

- (2) 12" (305mm) Heavy-Duty Stainless Steel File Drawers
- Slate Blue Drawer Pulls

Countertop Height: 31.5"H (800mm)

Overall Dimensions: 34.5"H (876mm) x 73"W (1855mm) x 28.6"D (727mm)



**SXRK32MW6**

**SXRK32MW6** includes:

73" (1855mm) Locking Mobile WorkCenter with 27" (686mm) Knee-hole & 3" (76mm) Total Lock Casters  
 Laminate Countertop with 3" (76mm) Backsplash (Wilsonart Carbon EV)  
 Pencil Drawer & Knee-hole Connecting Brackets

Left Mobile Bay (Keylock/2 keys)

- (2) 3" (76mm) Polymer Drawers
- (1) 6" (152mm) Polymer Drawer
- (1) 9" (229mm) Polymer Drawer

Right Mobile Bay (Keylock/2 keys)

- (2) 3" (76mm) Polymer Drawers
- (1) 6" (152mm) Polymer Drawer
- (1) 9" (229mm) Polymer Drawer

Slate Blue Drawer Pulls

Countertop Height: 31.5"H (800mm)

Overall Dimensions: 34.5"H (876mm) x 73"W (1855mm) x 28.6"D (727mm)



**SXRK41MW7**

**SXRK41MW7** includes:

73" Mobile Non-Locking Mobile WorkCenter 27" (686mm) Knee-hole & 3" (76mm) Total Lock Casters  
 Laminate Countertop with 3" (76mm) Backsplash (Wilsonart Carbon EV)  
 Pencil Drawer & Knee-hole Connecting Brackets

Left Mobile Bay

- (1) 3" (76mm) Polymer Drawers
- (2) 6" (152mm) Polymer Drawers
- (2) 9" (229mm) Polymer Drawers

Right Mobile Bay


- (3) Full-Extension Chrome Wire Shelves
- Right-Hinged 225° Clear Non-Locking Door

Slate Blue Drawer & Door Pulls


Countertop Height: 40.5"H (1029mm)

Overall Dimensions: 43.5"H (1105mm) x 73"W (1855mm) x 28.6"D (727mm)


**Standard Drawer Tote Specs**




**3" Drawer Tote**



**3" Drawer Tote**



**6" Drawer Tote**



**9" Drawer Tote**



**Weight Capacities:**  
 3" Tote = 25 lbs. (11.3Kg)  
 6" & 9" Totes = 40 lbs. (18.1 Kg)

**Heavy Duty Drawer Specs**



Above: Double-Wide 3" Heavy Duty Drawer

Left: Single-Wide Heavy Duty Drawers

Inside Dimensions Metro Part #	Inches			Millimeters			Wt Capacity	
	H	W	D	H	W	D	Lbs	Kg
<b>Single-Wide HD Drawers</b>								
3" Heavy Duty Dwr	2.3	16.3	19.2	58	414	487	100	45.4
6" Heavy Duty Dwr	5.3	16.3	19.2	134	414	487	100	45.4
9" Heavy Duty Dwr	8.3	16.3	19.2	210	414	487	100	45.4
12" Heavy Duty Dwr	11.5	16.3	19.2	292	414	487	100	45.4
<b>Double-Wide HD Drawers</b>								
3" Heavy Duty Dwr	2.3	36.3	19.2	58	922	487	150	68
6" Heavy Duty Dwr	5.3	36.3	19.2	134	922	487	150	68
9" Heavy Duty Dwr	8.3	36.3	19.2	210	922	487	150	68
12" Heavy Duty Dwr	11.5	36.3	19.2	292	922	487	150	68

**Construction:** The drawer interior is composed of 16 and 18 gauge solid stainless steel. The drawer slides are full-extension, zinc-plated steel with chrome-plated ball bearings, and attach to qwikSLOT inserts. They have side lever releases and detents in the closed position. The slide travel distance is 20.75". The drawer front surface is composed of molded ABS polymer, colored in Metro Taupe. Snap-on handles ("pulls") attach above the top recessed area and are available in a variety of colors. Optional self-adhesive label holders positioned in the recessed area under the pulls are also offered.



**SXRMS33EBS1**

**SXRMS33EBS1** includes:

Single-Wide Non-Locking Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
 (2) 12" (305mm) Heavy-Duty Stainless Steel File Drawers  
 Taupe Drawer Pulls  
Countertop Height: 33.5"H (1003mm)  
Overall Dimensions: 36.5"H (927mm) x 23.6"W (600mm) x 28.6"D (727mm)



**SXRMS36EBS2**

**SXRMS36EBS2** includes:

Single-Wide Non-Locking Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
 (2) 3" (76mm) Polymer Drawers  
 (2) 6" (152mm) Polymer Drawers  
 (1) 9" (229mm) Polymer Drawer  
 Taupe Drawer Pulls  
Countertop Height: 36.5"H (927mm)  
Overall Dimensions: 39.5"H (1003mm) x 23.6"W (600mm) x 28.6"D (727mm)



**SXRMS36EBS3**

**SXRMS36EBS3** includes:

Single-Wide Locking Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
 (3) Polymer Shelves  
 Left-Hinged 225° Solid Locking Door (Keylock/2 keys)  
 Taupe Door Pull  
Countertop Height: 36.5"H (927mm)  
Overall Dimensions: 39.5"H (1003mm) x 22.7"W (577mm) x 28.6"D (727mm)



**SXRMSD33EBS1**

**SXRMSD33EBS1** includes:

Double-Wide Non-Locking Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
 (4) 12" (305mm) Heavy-Duty Stainless Steel File Drawers  
 Taupe Drawer Pulls  
Countertop Height: 33.5"H (1003mm)  
Overall Dimensions: 36.5"H (927mm) x 43.4"W (1103mm) x 28.6"D (727mm)



**SXR MWD36EBS2**

**SXR MWD36EBS2** includes:

Double-Wide Non-Locking Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
 (4) 3" (76mm) Polymer Drawers  
 (4) 6" (152mm) Polymer Drawers  
 (2) 9" (229mm) Polymer Drawer  
 Taupe Drawer Pulls

Countertop Height: 36.5"H (927mm)

Overall Dimensions: 39.5"H (1003mm) x 43.4"W (1103mm) x 28.6"D (727mm)



**SXR MWD36EBS3**

**SXR MWD36EBS3** includes:

Double-Wide Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
 (6) Polymer Shelves  
 Left-Hinged 225° Solid Locking Door (Keylock/2 keys)  
 Right-Hinged 225° Solid Locking Door (Keylock/2 keys)  
 Taupe Door Pulls

Countertop Height: 36.5"H (927mm)

Overall Dimensions: 39.5"H (1003mm) x 43.4"W (1103mm) x 28.6"D (727mm)



**SXR MWT36EBS4**

**SXR MWT36EBS4** includes:

Triple-Wide Mobile WorkCenter Shell with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
Left Mobile Bay  
 (3) Single-Wide Full-Extension Chrome Wire Shelves  
 Left-Hinged 225° Solid Locking Door (Keylock/2 keys)

Right Mobile Bay

(3) Double-Wide Full-Extension Chrome Wire Shelves  
 Center-Closing 225° Solid Locking Door Set (Keylock/2 keys)  
 Taupe Door Pulls

Countertop Height: 36.5"H (927mm)

Overall Dimensions: 39.5"H (1003mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXR MWT36EBS6**

**SXR MWT36EBS6** includes:

Triple-Wide Mobile WorkCenter with 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

(2) 3" (76mm) Polymer Drawers  
 (2) 6" (152mm) Polymer Drawers  
 (1) 9" (229mm) Polymer Drawer

Right Mobile Bay

(2) 3" (76mm) Double-Wide Heavy Duty Stainless Steel Drawers  
 (2) 6" (152mm) Double-Wide Heavy Duty Stainless Steel Drawers  
 (1) 9" (229mm) Double-Wide Heavy Duty Stainless Steel Drawer  
 Taupe Drawer Pulls

Countertop Height: 36.5"H (927mm)

Overall Dimensions: 39.5"H (1003mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXRMTWK36EBS2**

**SXRMTWK36EBS2** includes:

Triple-Wide Non-Locking Mobile WorkCenter with 35" (890mm) Right-Side Knee-hole, Pencil Drawer & 5" (127mm) Diameter Polymer Casters (3-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

- (2) 3" (76mm) Polymer Drawers
- (2) 6" (152mm) Polymer Drawers
- (1) 9" (229mm) Polymer Drawer
- Taupe Drawer Pulls

Countertop Height: 36.5"H (927mm)

Overall Dimensions: 39.5"H (1003mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXRMTWK36EBS3**

**SXRMTWK36EBS3** includes:

Triple-Wide Mobile WorkCenter with 35" (890mm) Right-Side Knee-hole, Pencil Drawer & 5" (127mm) Diameter Polymer Casters (3-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

- (3) Polymer Shelves
- Left-Hinged 225° Solid Locking Door (Keylock/2 keys)
- Taupe Door & Drawer Pulls

Countertop Height: 36.5"H (927mm)

Overall Dimensions: 39.5"H (1003mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXRMTWK39EBS2**

**SXRMTWK39EBS2** includes:

Triple-Wide Non-Locking Mobile WorkCenter with 35" (890mm) Right-Side Knee-hole, Pencil Drawer & 5" (127mm) Diameter Polymer Casters (3-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

- (3) 3" (76mm) Polymer Drawers
- (2) 6" (152mm) Polymer Drawers
- (1) 9" (229mm) Polymer Drawer
- Taupe Drawer Pulls

Countertop Height: 39.5"H (1004mm)

Overall Dimensions: 42.5"H (1080mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXRMTWK39EBS3**

**SXRMTWK39EBS3** includes:

Triple-Wide Mobile WorkCenter with 35" (890mm) Right-Side Knee-hole, Pencil Drawer & 5" (127mm) Diameter Polymer Casters (3-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

- (3) Polymer Shelves
- Left-Hinged 225° Solid Locking Door (Keylock/2 keys)
- Taupe Door & Drawer Pulls

Countertop Height: 39.5"H (1004mm)

Overall Dimensions: 42.5"H (1080mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXRMSK33EBS1**

**SXRMSK33EBS1** includes:

73" (1855mm) Non-Locking Mobile WorkCenter with 27" (686mm) Center Knee-hole & 5" (127mm) Diameter Polymer Casters (4-Toe Brake, 4-Swivel)  
1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
Pencil Drawer & Knee-hole Connecting Brackets

Left Mobile Bay

(2) 12" (305mm) Heavy-Duty Stainless Steel File Drawers

Right Mobile Bay

(2) 12" (305mm) Heavy-Duty Stainless Steel File Drawers

Taupe Drawer Pulls

Countertop Height: 33.5"H (1003mm)

Overall Dimensions: 36.5"H (927mm) x 73"W (1855mm) x 28.6"D (727mm)



**SXRMSK39EBS2**

**SXRMSK39EBS2** includes:

73" (1855mm) Non-Locking Mobile WorkCenter with 27" (686mm) Center Knee-hole & 5" (127mm) Diameter Polymer Casters (4-Toe Brake, 4-Swivel)  
1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
Pencil Drawer & Knee-hole Connecting Brackets

Left Mobile Bay

(3) 3" (76mm) Polymer Drawers

(2) 6" (152mm) Polymer Drawers

(1) 9" (229mm) Polymer Drawer

Right Mobile Bay

(3) 3" (76mm) Polymer Drawers

(2) 6" (152mm) Polymer Drawers

(1) 9" (229mm) Polymer Drawer

Taupe Drawer Pulls

Countertop Height: 39.5"H (1004mm)

Overall Dimensions: 42.5"H (1080mm) x 73"W (1855mm) x 28.6"D (727mm)



**SXRMSK39EBS3**

**SXRMSK39EBS3** includes:

73" (1855mm) Mobile WorkCenter with 27" (686mm) Center Knee-hole & 5" (127mm) Diameter Polymer Casters (4-Toe Brake, 4-Swivel)  
1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash  
Pencil Drawer & Knee-hole Connecting Brackets

Left Mobile Bay

(3) Polymer Shelves

Left-Hinged 225° Solid Locking Door (Keylock/2 keys)

Right Mobile Bay

(3) 3" Polymer Drawers

Right-Hinged 225° Solid Locking Door (Keylock/2 keys)

Taupe Drawer & Door Pulls

Countertop Height: 39.5"H (1004mm)

Overall Dimensions: 42.5"H (1080mm) x 73"W (1855mm) x 28.6"D (727mm)



**SXRMS36EOH2**

**SXRMS36EOH2** includes:

Single-Wide Non-Locking Mobile WorkCenter Shell with Overhead Cabinet & 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Mobile Bay

(2) 3" (76mm) Polymer Drawers  
 (2) 6" (152mm) Polymer Drawers  
 (1) 9" (229mm) Polymer Drawer

Overhead Cabinet (27"H (686mm) Interior)

(2) Polymer Shelves  
 Left-Hinged 225° Clear Locking Door (Keylock/2 keys)  
 Taupe Drawer & Door Pulls

Overhead Uprights & Crossbars

Countertop Height: 36.5"H (927mm)

Countertop to Overhead Clearance: 16" (406mm)

Overall Dimensions: 81.5"H (2070mm) x 23.6"W (600mm) x 28.6"D (727mm)



**SXRMS36EOH3**

**SXRMS36EOH3** includes:

Single-Wide Mobile WorkCenter Shell with Overhead Cabinet & 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Mobile Bay

(3) Polymer Shelves  
 Left-Hinged 225° Solid Locking Door (Keylock/2 keys)  
Overhead Cabinet (27"H (686mm) Interior)

(2) Polymer Shelves

Left-Hinged 225° Clear Locking Door (Keylock/2 keys)

Taupe Door Pulls

Overhead Uprights & Crossbars

Countertop Height: 36.5"H (927mm)

Countertop to Overhead Clearance: 16" (406mm)

Overall Dimensions: 81.5"H (2070mm) x 23.6"W (600mm) x 28.6"D (727mm)



**SXRMSD36EOH2**

**SXRMSD36EOH2** includes:

Double-Wide Mobile WorkCenter Shell with Overhead Cabinet & 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)  
 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

(2) 3" (76mm) Polymer Drawers  
 (2) 6" (152mm) Polymer Drawers  
 (1) 9" (229mm) Polymer Drawer

Right Mobile Bay

(2) 3" (76mm) Polymer Drawers  
 (2) 6" (152mm) Polymer Drawers  
 (1) 9" (229mm) Polymer Drawer

Double-Wide Overhead Cabinet (27"H (686mm) Interior)

(2) Double-Wide Polymer Shelves

Center-Closing 225° Clear Locking Door Set (Keylock/2 keys)

Taupe Drawer & Door Pulls

Overhead Uprights & Crossbars

Countertop Height: 36.5"H (927mm)

Countertop to Overhead Clearance: 16" (406mm)

Overall Dimensions: 81.5"H (2070mm) x 43.4"W (1103mm) x 28.6"D (727mm)



**SXR MWD36EOH3**

**SXR MWD36EOH3** includes:

Double-Wide Mobile WorkCenter Shell with Overhead Cabinet & 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

(3) Polymer Shelves

Left-Hinged 225° Solid Locking Door (Keylock/2 keys)

Right Mobile Bay

(3) 3" Polymer Drawers

Right-Hinged 225° Solid Locking Door (Keylock/2 keys)

Double-Wide Overhead Cabinet (27" (686mm) Interior)

(2) Double-Wide Polymer Shelves

Center-Closing 225° Clear Locking Door Set (Keylock/2 keys)

Taupe Door Pulls

Overhead Uprights & Crossbars

Countertop Height: 36.5"H (927mm)

Countertop to Overhead Clearance: 16" (406mm)

Overall Dimensions: 81.5"H (2070mm) x 43.4"W (1103mm) x 28.6"D (727mm)



**SXR MWT36EOH6**

**SXR MWT36EOH6** includes:

Triple-Wide Non-Locking Mobile WorkCenter Shell with Overhead Cabinet & 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

(2) 3" (76mm) Polymer Drawers

(2) 6" (152mm) Polymer Drawers

(1) 9" (229mm) Polymer Drawer

Right Mobile Bay

(2) 3" (76mm) Double-Wide Heavy Duty Stainless Steel Drawers

(2) 6" (152mm) Double-Wide Heavy Duty Stainless Steel Drawers

(1) 9" (229mm) Double-Wide Heavy Duty Stainless Steel Drawer

Triple-Wide Overhead Cabinet (27" (686mm) Interior)

Left Bay

(2) Single-Wide Polymer Shelves

Left-Hinged 225° Clear Locking Door (Keylock/2 keys)

Right Bay

(2) Double-Wide Polymer Shelves

Center-Closing 225° Clear Locking Door Set (Keylock/2 keys)

Taupe Drawer & Door Pulls

Overhead Uprights & Crossbars

Countertop Height: 36.5"H (927mm)

Countertop to Overhead Clearance: 16" (406mm)

Overall Dimensions: 82"H (2070mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXRMT36EOH3**

**SXRMT36EOH3** includes:

Triple-Wide Mobile WorkCenter Shell with Overhead Cabinet & 5" (127mm) Diameter Polymer Casters (2-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

(3) Polymer Shelves

Left-Hinged 225° Solid Locking Door (Keylock/2 keys)

Center Mobile Bay

(3) Polymer Shelves

Left-Hinged 225° Solid Locking Door (Keylock/2 keys)

Right Mobile Bay

(3) Polymer Shelves

Right-Hinged 225° Solid Locking Door (Keylock/2 keys)

Triple-Wide Overhead Cabinet (27" (686mm) Interior)

Left Bay

(2) Single-Wide Polymer Shelves

Left-Hinged 225° Clear Locking Door (Keylock/2 keys)

Right Bay

(2) Double-Wide Polymer Shelves

Center-Closing 225° Clear Locking Door Set (Keylock/2 keys)

Taupe Door Pulls

Overhead Uprights & Crossbars

Countertop Height: 36.5"H (927mm)

Countertop to Overhead Clearance: 16" (406mm)

Overall Dimensions: 81.5"H (2070mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXRMTK36EOH2**

**SXRMTK36EOH2** includes:

Triple-Wide Non-Locking Mobile WorkCenter Shell with 35" (890mm) Right-Side Kneehole, Pencil Drawer & 5" (127mm) Diameter Polymer Casters (3-Toe Brake, 2-Swivel)

1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

(2) 3" Polymer Drawers

(2) 6" Polymer Drawers

(1) 9" Polymer Drawer

Triple-Wide Overhead Cabinet (27" (686mm) Interior)

Left Bay

(2) Single-Wide Polymer Shelves

Left-Hinged 225° Clear Locking Door (Keylock/2 keys)

Right Bay

(2) Double-Wide Polymer Shelves

Center-Closing 225° Clear Locking Door Set (Keylock/2 keys)

Taupe Drawer & Door Pulls

Overhead Uprights & Crossbars

Countertop Height: 36.5"H (927mm)

Countertop to Overhead Clearance: 16" (406mm)

Overall Dimensions: 81.5"H (2070mm) x 63.2"W (1606mm) x 28.6"D (727mm)



**SXRMWTK36EOH3**

**SXRMWTK36EOH3** includes:

Triple-Wide Mobile WorkCenter Shell with 35" (890mm) Right-Side Knee-hole, Pencil Drawer & 5" (127mm) Diameter Polymer Casters (3-Toe Brake, 2-Swivel) 1" Thick Black Epoxy Countertop with 3" (76mm) Backsplash

Left Mobile Bay

(3) Polymer Shelves

Left-Hinged 225° Solid Locking Door (Keylock/2 keys)

Triple-Wide Overhead Cabinet (27" (686mm) Interior)

Left Bay

(2) Single-Wide Polymer Shelves

Left-Hinged 225° Clear Locking Door (Keylock/2 keys)

Right Bay

(2) Double-Wide Polymer Shelves

Center-Closing 225° Clear Locking Door Set (Keylock/2 keys)

Taupe Door Pulls

Overhead Uprights & Crossbars

Countertop Height: 36.5"H (927mm)

Countertop to Overhead Clearance: 16" (406mm)

Overall Dimensions: 81.5"H (2070mm) x 63.2"W (1606mm) x 28.6"D (727mm)

***Undercounter Unit (Phenolic Top)***



**SXRMWS29TS2**

**SXRMWS29TS2** includes:

Single-Wide Non-Locking Mobile Workcenter Shell with 3" (76mm)

Diameter Total-Lock Casters

1" Thick Black Phenolic Island Countertop (no backsplash)

(1) 4.5" (114mm) Polymer Drawer

(1) 3" (76mm) Polymer Drawer





(1) 7.5" (191mm) Polymer Drawer

(1) 6" (152mm) Polymer Drawer

Black Drawer Pulls

Overall Dimensions: 28.8"H (732mm) x 22.7"W (577mm) x 24.9"D (633mm)

**Standard Drawer Tote Specs**



**Weight Capacities:**  
 3" Tote = 25 lbs. (11.3Kg)  
 6" & 9" Totes = 40 lbs. (18.1 Kg)

**Heavy Duty Drawer Specs**



**Above:** Double-Wide 3" Heavy Duty Drawer

**Left:** Single-Wide Heavy Duty Drawers

Inside Dimensions Metro Part #	Inches			Millimeters			Wt Capacity	
	H	W	D	H	W	D	Lbs	Kg
<b>Single-Wide HD Drawers</b>								
3" Heavy Duty Dwr	2.3	16.3	19.2	58	414	487	100	45.4
6" Heavy Duty Dwr	5.3	16.3	19.2	134	414	487	100	45.4
9" Heavy Duty Dwr	8.3	16.3	19.2	210	414	487	100	45.4
12" Heavy Duty Dwr	11.5	16.3	19.2	292	414	487	100	45.4
<b>Double-Wide HD Drawers</b>								
3" Heavy Duty Dwr	2.3	36.3	19.2	58	922	487	150	68
6" Heavy Duty Dwr	5.3	36.3	19.2	134	922	487	150	68
9" Heavy Duty Dwr	8.3	36.3	19.2	210	922	487	150	68
12" Heavy Duty Dwr	11.5	36.3	19.2	292	922	487	150	68

**Construction:** The drawer interior is composed of 16 and 18 gauge solid stainless steel. The drawer slides are full-extension, zinc-plated steel with chrome-plated ball bearings, and attach to qwikSLOT inserts. They have side lever releases and detents in the closed position. The slide travel distance is 20.75". The drawer front surface is composed of molded ABS polymer, colored in Metro Taupe. Snap-on handles ("pulls") attach above the top recessed area and are available in a variety of colors. Optional self-adhesive label holders positioned in the recessed area under the pulls are also offered.

## STARSYS™ MOBILE, MODULAR WORKCENTERS

### *Modular System . . . Flexible Space*

Metro Starsys™ flexible lab furniture is a polymer-based, highly-adaptive system, that enhances lab productivity and space efficiency by empowering the end user with the ability to easily reconfigure and move units around the lab.

The Starsys Modular System creates flexible space. No matter what you need — open, closed, short, tall, mobile or stationary, preconfigured workcenters — Starsys has the solution for you.

With Starsys, you choose from an array of modular units to fit the space and work specifications. It's that easy. And that flexible.

### Benefits:

- Facilitates fast, efficient product or furniture installation.
- The ability to easily reconfigure lab furniture layout as needs change.
- The ability to easily reconfigure unit design.
- Durable, easy-to-clean, polymer enclosures that won't chip, peel, dent, or rust — ever.
- Rigid anodized aluminum and metal substructure, combined with advanced polymers, provides a robust and stable work environment for even the most sensitive equipment.
- Promote cGMP Protocol — Microban®\* antimicrobial product protection is built into every Starsys Unit inhibiting the growth of mold, mildew and bacteria.
- Advanced Security Options available — including electronic locks, proximity card readers, as well as mechanical keylocks.

\*MICROBAN® and the MICROBAN® symbol are registered trademarks of the Microban Products Company, Huntersville, NC.



**Double-wide  
Mobile WorkCenter  
with Reagent Shelving**



**Double-wide  
Mobile WorkCenter with  
Epoxy Top, Sink and  
Overhead Cabinets**



**Triple-wide  
Mobile WorkCenter  
with Kneewell,  
Overhead Storage  
and Stainless  
Steel Top**

To configure a Starsys unit to your exact specifications, please visit [www.metroconfigurator.com](http://www.metroconfigurator.com)



**InterMetro Industries Corporation**  
North Washington Street  
Wilkes-Barre, PA 18705  
[www.metro.com](http://www.metro.com)



## Specifications

Starsys™ epoxy resin worksurfaces are generally 1" (25mm) thick. Unless otherwise specified, worksurfaces are monolithic and molded from a modified epoxy resin, with a uniform mixture throughout the top forming a thermosetting material. Back or side splashes are shipped separately and affixed on site to epoxy top. Worksurfaces have a smooth, non-glare finish.

Phenolic countertops are 1" (25mm) thick with back and side splash options. The laboratory-grade tops are made of an epoxy resin composite (thermoset resins with cellulose fiber reinforcement), are non-porous, highly resistant to a wide variety of aggressive chemicals, easy to clean, disinfect and maintain.

Stainless steel countertops are 1.25" (32mm) thick with Type 304 series 18-gauge fabricated stainless steel with homopolymer polypropylene cores. Available with backsplashes or "island style."

Refer to Starsys chemical resistance report for a complete description of Starsys material specifications. A copy of the report can be found at [www.metro.com](http://www.metro.com) or by contacting your local Metro representative.

## Drawers

### Standard Drawers

Description	Inside Dimensions Height/Width/Length		Weight Capacity		Cat. No.
	(in.)	(mm)	(lbs.)	(kg)	
3" (76mm)	2 <sup>3</sup> / <sub>8</sub> x16 <sup>3</sup> / <sub>8</sub> x15 <sup>7</sup> / <sub>8</sub>	60x416x403	25	11.4	SXRS3
4.5" (114mm)	3 <sup>7</sup> / <sub>8</sub> x16 <sup>3</sup> / <sub>8</sub> x15 <sup>7</sup> / <sub>8</sub>	98x416x403	25	11.4	SXRS4.5
6" (152mm)	5 <sup>3</sup> / <sub>8</sub> x16 <sup>1</sup> / <sub>4</sub> x15 <sup>3</sup> / <sub>4</sub>	136x413x400	40	18.2	SXRS6
7.5" (191mm)	6 <sup>7</sup> / <sub>8</sub> x16 <sup>1</sup> / <sub>4</sub> x15 <sup>3</sup> / <sub>4</sub>	175x413x400	40	18.2	SXRS7.5
9" (230mm)	8 <sup>3</sup> / <sub>8</sub> x16 <sup>1</sup> / <sub>4</sub> x15 <sup>3</sup> / <sub>4</sub>	213x413x400	40	18.2	SXRS9
10.5" (267mm)	9 <sup>7</sup> / <sub>8</sub> x16 <sup>1</sup> / <sub>4</sub> x15 <sup>3</sup> / <sub>4</sub>	251x413x400	40	18.2	SXRS10.5
12" (305mm)	11 <sup>3</sup> / <sub>8</sub> x16 <sup>1</sup> / <sub>4</sub> x15 <sup>3</sup> / <sub>4</sub>	289x413x400	40	18.2	SXRS12



### Heavy-Duty Drawers

Description	Inside Dimensions Height/Width/Length		Mobile WorkCenters with Overhead Cabinets		Weight Capacity Stationary Applications		Mobile WorkCenters without Overhead Cabinets		Cat. No.
	(in.)	(mm)	(lbs.)	(kg)	(lbs.)	(kg)	(lbs.)	(kg)	
<b>Single-Wide Drawers</b>									
3" (76mm)	2 <sup>1</sup> / <sub>4</sub> x16 <sup>3</sup> / <sub>8</sub> x19 <sup>3</sup> / <sub>16</sub>	57x416x487	100	45.4	100	45.4	100	45.4	SXRS3HD
6" (152mm)	5 <sup>1</sup> / <sub>4</sub> x16 <sup>3</sup> / <sub>8</sub> x19 <sup>3</sup> / <sub>16</sub>	133x416x487	100	45.4	100	45.4	100	45.4	SXRS6HD
9" (230mm)	8 <sup>1</sup> / <sub>4</sub> x16 <sup>3</sup> / <sub>8</sub> x19 <sup>3</sup> / <sub>16</sub>	209x416x487	100	45.4	100	45.4	100	45.4	SXRS9HD
12" (305mm)*	11 <sup>1</sup> / <sub>4</sub> x16 <sup>3</sup> / <sub>8</sub> x19 <sup>3</sup> / <sub>16</sub>	285x416x487	100	45.4	100	45.4	100	45.4	SXRS12HD
<b>Double-Wide Drawers</b>									
3" (76mm)	2 <sup>1</sup> / <sub>8</sub> x36 <sup>1</sup> / <sub>4</sub> x19 <sup>3</sup> / <sub>16</sub>	54x920x487	150	68.0	150	68.0	150	68.0	SXRD3HD
6" (152mm)	5 <sup>1</sup> / <sub>8</sub> x36 <sup>1</sup> / <sub>4</sub> x19 <sup>3</sup> / <sub>16</sub>	130x920x487	150	68.0	150	68.0	150	68.0	SXRD6HD
9" (230mm)	8 <sup>1</sup> / <sub>8</sub> x36 <sup>1</sup> / <sub>4</sub> x19 <sup>3</sup> / <sub>16</sub>	206x920x487	150	68.0	150	68.0	150	68.0	SXRD9HD
12" (305mm)	11 <sup>1</sup> / <sub>8</sub> x36 <sup>1</sup> / <sub>4</sub> x19 <sup>3</sup> / <sub>16</sub>	283x920x487	150	68.0	150	68.0	150	68.0	SXRD12HD

\*Note: Single wide 12" (305mm) file drawers include hanging file accessory.

To specify color for drawer pull, order following catalog number for each drawer:

Description	Cat. No.
Black	SXRDP-BK
Taupe	SXRDP-TP
Slate Blue	SXRDP-BL
Green	SXRDP-GR
Code Blue	SXRDP-CB
Yellow	SXRDP-YL
Violet	SXRDP-VL
Orange	SXRDP-OR
Red	SXRDP-RE
White	SXRDP-WH
Pink	SXRDP-PK

Note: 2 drawer pulls are required for double-wide heavy-duty drawers.  
Note: Actual color images can be seen in the Commercial Solutions Guide or at [www.metro.com](http://www.metro.com).



Single bin and Storage level access modes available.

### Drawer/Door Label Holders

10 per package.  
3/4" x 11" (19 X 279mm)  
Cat. No. **SXRLABKIT**



Starsys Carts are built to exact customer specifications, therefore, they are non-cancelable, non-returnable.

## Mobile WorkCenter Drawer and Door Comparison

To select an appropriate drawer or door configuration for your application, consult your Metro representative or visit [www.metroconfigurator.com](http://www.metroconfigurator.com)

### Starsys™ Drawers

Starsys offers two types of drawers to address varying user requirements. Below is a feature comparison to help determine the appropriate drawer system for a desired application.

Feature	Polymer Drawers	Heavy-Duty Drawers
Light to medium weight capacity (25-40 lbs.)	X	
Heavy weight capacity (up to 150 lbs.)		X
Removable tote	X	
Inner panel compatible	X	
Slotted post insert compatible		X
Ball bearing slide		X
Full extension	X	X
Stainless steel interior		X
Polymer interior	X	
Easily removed/reconfigured	X	
Adjustable dividers	X	X
Lockable	X	X
Optional lock cover	X	
Label kits	X	X
Accommodate hanging files		X



Removable drawer totes allow for exchange or convenient access outside of cart



Heavy-Duty Drawer

### Starsys™ Doors

Starsys offers two types of doors, hinged and tambour. The hinged doors are available in a choice of solid or clear, locking or non-locking. Each offers unique benefits to address varying user requirements. Below is a feature comparison to help determine the appropriate door selection for a desired application.

Feature	Hinged Solid Door	Hinged Clear Door	Rollup Tambour Door	Benefit
Visual Security	X		X	Puts inventory out of plain site
Hides clutter	X		X	Enhances department aesthetics
Visual inventory		X		Protection is not compromised when looking for supplies
Lockable	X	X	X	Added security
Space saving			X	Does not encroach on work space
Door swing overlaps adjacent storage spaces	X	X		Promotes door closure to protect supplies
90 degree hinge	X	X		Prevents door from overlapping adjacent space
270 degree hinge	X	X		Allows door to swing to side of cabinet or cart
Polymer	X	X	X	Impact and corrosion protection
Left or right hinge	X	X		Provides appropriate swing direction
Center closing	X	X		Allows for double wide storage compartments
Door label holders	X	X		Provides clean replenishable labeling option
Non-locking	X	X		Lower cost alternative



Tall Clear Locking Center-Closing Door Set for units with double-wide wire shelves



Clear Locking Center-Closing Door Set on Overhead and solid single wide non-locking door on mobile base



Tall Clear Single Wide locking door on unit with poly shelves

To configure a Starsys unit to your exact specifications, please visit [www.metroconfigurator.com](http://www.metroconfigurator.com)

Mobile WorkCenter Overhead Cabinets\*  
Cabinets include inner panels and mounting rails.



**Double-wide Mobile WorkCenter with Epoxy Top, Sink and Overhead Cabinets**

Description	Height/Depth/Width Overall Dimensions		Inside Dimensions		Cat. No.
	(in.)	(mm)	(in.)	(mm)	
<b>Single-Wide Cabinets</b>					
24" (610mm)	26 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	670x343x546	24 <sup>5</sup> / <sub>8</sub> x12x17	625x341x432	SXRSOH24P
27" (686mm)	29 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	746x343x546	27 <sup>5</sup> / <sub>8</sub> x12x17	702x341x432	SXRSOH27P
30" (760mm)	32 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	822x343x546	30 <sup>5</sup> / <sub>8</sub> x12x17	778x341x432	SXRSOH30P
33" (838mm)	35 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	898x343x546	33 <sup>5</sup> / <sub>8</sub> x12x17	854x341x432	SXRSOH33P
36" (914mm)	38 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	973x343x546	36 <sup>5</sup> / <sub>8</sub> x12x17	930x341x432	SXRSOH36P
<b>Double-Wide Cabinets</b>					
24" (610mm)	26 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	670x343x1051	24 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	625x254x933	SXRDOH24P
27" (686mm)	29 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	746x343x1051	27 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	702x254x933	SXRDOH27P
30" (760mm)	32 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	822x343x1051	30 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	778x254x933	SXRDOH30P
33" (838mm)	35 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	898x343x1051	33 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	854x254x933	SXRDOH33P
36" (914mm)	38 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	973x343x1051	36 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	930x254x933	SXRDOH36P
<b>Shelves for Single-Wide Cabinets †</b>					
Wire Shelf	1 <sup>3</sup> / <sub>8</sub> x11x17 <sup>15</sup> / <sub>16</sub>	35x279x456			SXRSOHWE
Poly Shelf	1 <sup>1</sup> / <sub>2</sub> x12x18	38x305x457			SXRSOHPS
<b>Shelves for Double-Wide Cabinets ††</b>					
Wire Shelf	1 <sup>3</sup> / <sub>8</sub> x11x36 <sup>11</sup> / <sub>16</sub>	35x279x932			SXRDOHWE
Poly Shelf	1 <sup>1</sup> / <sub>2</sub> x12x37 <sup>5</sup> / <sub>8</sub>	38x305x956			SXRDOHPS

† Load rating for single-wide shelves, both wire and poly, is 60 lbs. (27kg).  
 †† Load rating for double-wide shelves, both wire and poly, is 120 lbs. (54kg).  
 \* Triple-wide overhead cabinets are only offered as part of Mobile WorkCenter offering. Contact your Metro representative for more information.

Overhead Accessories

Description	Height/Depth/Width		Cat. No.
	(in.)	(mm)	
Light for Single-Wide Overhead Cabinets	1.50x5.82x19.75	38x148x502	SXRSOHL
Light for Double-Wide Overhead Cabinets	1.50x5.82x39.55	38x148x1005	SXRDOHL
Light for Triple-Wide Overhead Cabinets	1.50x5.82x59.35	38x148x1507	SXRTOHL
Overhead Light (No Valance)	1.32x4.63x14.25	34x118x362	SXROHLT
Overhead Back Filler Strip (1 Pair)*	2.875x.125x49.75	73x3x1264	SXROHFLR-49

\*One kit recommended for a single unit or for each end on a run of multiple units.



SXRLABKIT

Drawer/Door Label Holders

Description	Dimensions		Cat. No.
	Height	Length	
Label Kits (10 Pieces Per Kit)	3/4"	11" (19 x 279mm)	SXRLABKIT



**Single-Wide Overhead Cabinet with Solid Door**



**Double-Wide Overhead Cabinet with Wire Shelves (2) and Locking, Clear Doors**



**Double-Wide Overhead Cabinet with Polymer Shelves and Locking, Clear Doors**

Starsys Carts are built to exact customer specifications, therefore, they are non-cancelable, non-returnable.

## Mobile WorkCenter Overhead Cabinet and Reagent Shelving Support

Description	Cat. No.	
<b>Lateral Support for Overhead Cabinets — Select one per WorkCenter</b>		
Overhead Structure Single Wide Mobile WorkCenter		SXRSMWOHB
Overhead Structure Double Wide Mobile WorkCenter		SXRDMWOHB
Overhead Structure Triple Wide Mobile WorkCenter		SXRTMWOHB
Overhead Structure Single-Kneewell Mobile WorkCenter		SXRKMWOHB
<b>Vertical Supports for Overhead Cabinets — Select one pair per WorkCenter</b>		
24" (610mm) Base Interior with 24" (610mm) Overhead Interior	(1 Pair)	SXRMWUP24-24
24" (610mm) Base Interior with 27" (686mm) Overhead Interior	(1 Pair)	SXRMWUP24-27
24" (610mm) Base Interior with 30" (762mm) Overhead Interior	(1 Pair)	SXRMWUP24-30
24" (610mm) Base Interior with 33" (838mm) Overhead Interior	(1 Pair)	SXRMWUP24-33
27" (686mm) Base Interior with 24" (610mm) Overhead Interior	(1 Pair)	SXRMWUP27-24
27" (686mm) Base Interior with 27" (686mm) Overhead Interior	(1 Pair)	SXRMWUP27-27
27" (686mm) Base Interior with 30" (762mm) Overhead Interior	(1 Pair)	SXRMWUP27-30
30" (762mm) Base Interior with 24" (610mm) Overhead Interior	(1 Pair)	SXRMWUP30-24
30" (762mm) Base Interior with 27" (686mm) Overhead Interior	(1 Pair)	SXRMWUP30-27
33" (838mm) Base Interior with 24" (610mm) Overhead Interior	(1 Pair)	SXRMWUP33-24
<b>Overhead Cabinet and Shelving Combination — Vertical Supports — Select one pair per WorkCenter</b>		
24" (610mm) Base 24" (610mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP24Q-24
24" (610mm) Base 27" (686mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP24Q-27
24" (610mm) Base 30" (762mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP24Q-30
24" (610mm) Base 33" (838mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP24Q-33
27" (686mm) Base 24" (610mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP27Q-24
27" (686mm) Base 27" (686mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP27Q-27
27" (686mm) Base 30" (762mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP27Q-30
30" (762mm) Base 24" (610mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP30Q-24
30" (762mm) Base 27" (686mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP30Q-27
33" (838mm) Base 24" (610mm) Overhead Interior & Slotted Inserts	(1 Pair)	SXRMWUP33Q-24

Contact the Metro Starsys team for dimension and packaged weight information.



**Mobile WorkCenter with Reagent Shelving**

\*Base unit working height is equal to Base interior height +9" [229mm] (based on 5" [127mm] casters)

Overhead exterior height is equal to Overhead interior height +2.5" [64mm]

## Mobile WorkCenter Reagent Shelving and Support

Description	Cat. No.	
<b>Reagent Shelving Options</b>		
Shelf Assembly — Single width		SXRSMWSHF
Shelf Assembly — Double width		SXRDMWSHF
Shelf Assembly — Triple width		SXRTMWSHF
Shelf Assembly — Single — Kneewell — Single		SXRKMWSHF
<b>Lateral Support for Reagent Shelving — Select one pair per WorkCenter</b>		
Cross Member — Single width	(1 Pair)	SXRSMWSB
Cross Member — Double width	(1 Pair)	SXRDMWSB
Cross Member — Triple width	(1 Pair)	SXRTMWSB
Cross Member — Single — Kneewell — Single	(1 Pair)	SXRKMWSB
<b>Vertical Supports with Slotted Inserts for Reagent Shelving — Select one per WorkCenter</b>		
24" (610mm) Base Interior with Slotted Inserts	(1 Pair)	SXRMWUP24Q
27" (686mm) Base Interior with Slotted Inserts	(1 Pair)	SXRMWUP27Q
30" (762mm) Base Interior with Slotted Inserts	(1 Pair)	SXRMWUP30Q
33" (838mm) Base Interior with Slotted Inserts	(1 Pair)	SXRMWUP33Q

Contact the Metro Starsys team for dimension and packaged weight information.

## Casters Options

All Mobile WorkCenters come standard with (2) Polymer Total-Lock front side casters and (2) Polymer Swivel back side casters.

Description	Cat. No.
(2) Front Polymer 5" (127mm) Total-Lock Casters	B5P-TL
(2) Rear Polymer 5" (127mm) Swivel Casters	B5P-SW

Additional caster options are available through the Metro Starsys Configurator or by contacting your Metro representative.

<p><b>3" (76mm) Total-Lock Caster:</b> Utilize 3" total-lock plate casters for limited mobile applications or when additional storage area is required.</p>  <p>B3P-TL</p>	<p><b>5" (127mm) Stainless Steel Caster:</b> Recommended for corrosive environments.</p>  <p>B5PBGSA</p>	<p><b>5" (127mm) Total-Lock, All-Polymer Caster:</b> Recommended for corrosive environments. All-polymer total-lock plate caster is an economical alternative to stainless casters.</p>  <p>B5PC-TL</p>	<p><b>Stabilizer/Leveling Caster:</b> Recommended for applications where benchtop equipment demands unit stability and/or a level worksurface. Engage caster foot for stability, release for mobility.</p>  <p>B3N-AS</p>
<p>5" (127mm) Swivel/Lock (Directional) Casters (available by request) Other specialty caster options are available, consult your Metro representative.</p>			

## Preconfigured Mobile WorkCenters with Black Epoxy Tops

All preconfigured units below come standard with taupe door/drawer pulls and polymer 5" (127mm) diameter casters (2-total lock, 2-swivel). Overhead units include 27" (686mm) cabinets with locking clear doors and two poly shelves. The countertops are 1" (25mm) thick black epoxy resin. Other top materials and accessories are available upon request. Contact your Metro representative or visit [www.metroconfigurator.com](http://www.metroconfigurator.com).

Description	Cat. No.	Overall Dimensions Width/Length/Height	Cat. No. With Overhead	Overall Dimensions Width/Length/Height
<b>Single-Wide Units</b>				
Single-Wide Heavy-Duty 12" (305mm) (2)	SXRMWS33EBS1	28.6x23.6x33.5		
Single-Wide with Standard Drawers	SXRMWS36EBS2	28.6x23.6x36.5	SXRMWS36EOH2	28.6x23.6x81.5
Single-Wide with Locking Doors and Poly Shelves	SXRMWS36EBS3	28.6x23.6x36.5	SXRMWS36EOH3	28.6x23.6x81.5
<b>Double-Wide Units</b>				
Double-Wide Heavy-Duty 12" (305mm) (2)	SXRMWD33EBS1	28.6x43.4x33.5		
Double-Wide with Standard Drawers	SXRMWD36EBS2	28.6x43.4x36.5	SXRMWD36EOH2	28.6x43.4x81.5
Double-Wide with Locking Doors and Poly Shelves	SXRMWD36EBS3	28.6x43.4x36.5	SXRMWD36EOH3	28.6x43.4x81.5
<b>Triple-Wide Units</b>				
Triple-Wide with S & D HD Drawers	SXRMWT36EBS6	28.6x63.2x36.5	SXRMWT36EOH6	28.6x63.2x81.5
Triple-Wide with Locking Doors and Wire Shelves	SXRMWT36EBS4	28.6x63.2x36.5		
Triple-Wide with Locking Doors and Poly Shelves			SXRMWT36EOH3	28.6x63.2x36.5
Triple Kneewell with Standard Drawers	SXRMWTK36EBS2	28.6x63.2x36.5	SXRMWTK36EOH2	28.6x63.2x81.5
Triple Kneewell with Locking Doors and Poly Shelves	SXRMWTK36EBS3	28.6x63.2x36.5	SXRMWTK36EOH3	28.6x63.2x81.5
Triple Kneewell with Standard Drawers	SXRMWTK39EBS2	28.6x63.2x39.5		
Triple Kneewell with Locking Doors and Poly Shelves	SXRMWTK39EBS3	28.6x63.2x39.5		
Single Kneewell with HD Drawers	SXRMWSK33EBS1	28.6x73x33.5		
Single Kneewell with Standard Drawers	SXRMWSK39EBS2	28.6x73x39.5		
Single Kneewell with Locking Doors and Poly Shelves	SXRMWSK39EBS3	28.6x73x39.5		



SXRMWS33EBS1



SXRMWS36EBS2



SXRMWS36EBS3



SXRMWD33EBS1



SXRMWD36EBS2



SXRMWD36EBS3



SXRMWT36EBS4



SXRMWT36EBS6

Starsys Carts are built to exact customer specifications, therefore, they are non-cancelable, non-returnable.

Starsys Mobile Workcenter models shown below include epoxy countertops. These models are also available with stainless steel, phenolic, solid surface and laminate tops. Contact your Metro representative for more information or visit [www.metroconfigurator.com](http://www.metroconfigurator.com)



SXRMTWK36EBS2



SXRMTWK36EBS3



SXRMTWK39EBS2



SXRMTWK39EBS3



SXRMTWSK33EBS1



SXRMTWSK39EBS2



SXRMTWSK39EBS3



SXRMTWS36EOH2



SXRMTWS36EOH3



SXRMTWD36EOH2

To configure a Starsys unit to your exact specifications, please visit [www.metroconfigurator.com](http://www.metroconfigurator.com)

Job \_\_\_\_\_

LAB SOLUTIONS

**METRO**

Starsys Mobile Workcenter models shown below include epoxy countertops. These models are also available with stainless steel, phenolic, solid surface and laminate tops. Contact your Metro representative for more information or visit [www.metroconfigurator.com](http://www.metroconfigurator.com)



SXR MWD36EOH3



SXR MWT36EOH6



SXR MWT36EOH3



SXR MWTk36EOH2



SXR MWTk36EOH3

**Starsys Carts are built to exact customer specifications, therefore, they are non-cancelable, non-returnable.**

All Metro Catalog Sheets are available on our Web Site: [www.metro.com](http://www.metro.com)



**InterMetro Industries Corporation**  
North Washington Street, Wilkes-Barre, PA 18705  
Phone: 570-825-2741  
Fax: 570-825-2852

L03-252  
Printed in U.S.A. Rev. 1/11

Information and specifications are subject to change without notice. Please confirm at time of order.

Copyright © 2011 InterMetro Industries Corp.

For Product Information:  
U.S. and Canada: 1.800.992.1776

Latin America: 1.561.961.1389  
Europe: +31.76.587.7550

Asia/Pacific: +65.6567.8003  
Middle East/Africa: +971.4.811.8286

