



Item # \_\_\_\_\_

Job \_\_\_\_\_



## CENTRAL STERILE PRODUCTIVITY BENCH

This productivity bench provides a clean, obstruction-free surface for the preparation of surgical instrument trays and surgical packs and an optional accessory grid for tools and materials required.

### Features Include:

- Adjustable legs, footrest or shelf to help reduce worker fatigue.
- Adjustable work surface height is from 31"-37" (787-940mm) on 1" (25mm) increments.
- Type 304 stainless steel work surface.
- Rounded corners to reduce injury and torn packs.
- Posts are HD 1 1/2" (38mm) diameter 16 gauge powder coated steel.
- Wide range of accessories to customize by need.



### Basic models consist of:

- Two short front posts and two tall rear posts with accessory slots.
- Stainless steel work surface.
- Support frame(s): bottom three-sided and top tie-bar.
- Choice of solid bottom shelf or foot rest frame.

### Basic Model with Stainless Steel Work Surface\*

| Model No. | Depth x Width |            | Description                               | Approx. Pkd. Wt. |      |
|-----------|---------------|------------|---|------------------|------|
|           | (in.)         | (mm)       |   | (lbs.)           | (kg) |
| PBH3072FS | 30 x 72       | 762 x 1829 | 72" (1829mm) Bench with Flat Bottom Shelf | 257              | 116  |
| PBH3072HS | 30 x 72       | 762 x 1829 | 72" (1829mm) Bench with H-Frame Footrest  | 238              | 107  |
| PBH3096FS | 30 x 96       | 762 x 2438 | 96" (2438mm) Bench with Flat Bottom Shelf | 291              | 131  |
| PBH3096HS | 30 x 96       | 762 x 2438 | 96" (2438mm) Bench with H-Frame Footrest  | 272              | 123  |

\*Actual work surface depth is 30" (762mm).  
Work surface weight capacity is 1,000 lbs. (453kg) evenly distributed.

## Accessories

### Stainless Steel Work Surface Mounted

| Model No. | Width x Length x Height |                 | Description             | Approx. Pkd. Wt. |      |
|-----------|-------------------------|-----------------|-------------------------|------------------|------|
|           | (in.)                   | (mm)            |                         | (lbs.)           | (kg) |
| WTD21C    | 14 1/4 x 28 x 5 3/4     | 362 x 711 x 146 | Economy Polymer Drawer  | 7                | 3    |
| WTD51S*   | 24 x 25 3/4 x 7 1/2     | 610 x 654 x 191 | Stainless Deluxe Drawer | 31               | 14   |

\*WTD51S is stackable.



**InterMetro Industries Corporation**  
North Washington Street  
Wilkes-Barre, PA 18705  
www.metro.com





# CENTRAL STERILE PRODUCTIVITY BENCH

## Accessories

### Post Mounted

| Model No.               | Width x Length x Height  |                 | Description                        | Approx. Pkd. Wt. |      |
|-------------------------|--|-----------------|------------------------------------|------------------|------|
|                         | (in.)  | (mm)            |                                    | (lbs.)           | (kg) |
| PBA-1260CQ <sup>1</sup> | 12 x 60 x 1 <sup>1</sup> / <sub>4</sub>  | 305 x 1525 x 31 | (A) Adjustable Cantilever Shelf    | 13               | 5.9  |
| PBA-60PQ <sup>2</sup>   | 1 <sup>1</sup> / <sub>2</sub> x 57 <sup>3</sup> / <sub>4</sub> x 1 <sup>1</sup> / <sub>2</sub> | 35 x 1467 x 41  | (B) Power Outlet Strip             | 4                | 1.8  |
| PBA-LA60 <sup>3</sup>   | 22 x 55 x 3  | 305 x 1397 x 76 | (C) Light Assembly                 | 26               | 11.7 |
| PBA-BBCQ                | 1 <sup>1</sup> / <sub>2</sub> x 57 <sup>3</sup> / <sub>4</sub> x 3                             | 13 x 1467 x 76  | Bin Rail                           | 10               | 4.5  |
| PBA-GPCQ <sup>4</sup>   | 33 x 57 <sup>3</sup> / <sub>4</sub>  | 1467 x 838      | Full Height Grid with Clips        | 21               | 9.5  |
| WG1860BR <sup>4†</sup>  | 18 x 60  | 457 x 1524      | (D) Half Height Grid               | 14               | 6.3  |
| 9985QS <sup>†</sup>     |  |                 | (E) Grid Clips (package of 4)      | .3               | .1   |
| 9985QSL <sup>†</sup>    |  |                 | (E) Grid Clip Locks (package of 4) | .3               | .1   |

<sup>†</sup>Requires 1 - 9985QS and 1 - 9985QSL.

### Grid Mounted

| Model No.             | Width x Length x Height   |                 | Description             | Approx. Pkd. Wt. |      |
|-----------------------|---|-----------------|-------------------------|------------------|------|
|                       | (in.)   | (mm)            |                         | (lbs.)           | (kg) |
| PBA-GSD <sup>5</sup>  | 9 x 18 <sup>1</sup> / <sub>2</sub>                                | 228 x 470       | (F) Grid Shelf          | 4                | 2    |
| PBA-MS <sup>5</sup>   | 16 <sup>3</sup> / <sub>4</sub> x 16 <sup>3</sup> / <sub>4</sub>   | 425 x 425       | (G) Monitor Shelf       | 6                | 3    |
| PBA-PFH               | 12 <sup>1</sup> / <sub>2</sub> x 6 <sup>1</sup> / <sub>4</sub>    | 317 x 159       | Hanging File Holder     | 3                | 1    |
| PBA-CHD               | 2 x 9   | 50 x 228        | (H) Catalog/File Holder | 3                | 1    |
| PBA-1BH               | 11 <sup>3</sup> / <sub>4</sub> x 3                                | 298 x 76        | (I) Short Bin Holder    | 2                | 1    |
| PBA-2BH               | 22 <sup>1</sup> / <sub>4</sub> x 3                                | 565 x 76        | Long Bin Holder         | 3                | 1    |
| MB30230B <sup>6</sup> | 5 <sup>1</sup> / <sub>2</sub> x 11 x 5                            | 140 x 279 x 127 | (J) Small Bin           | 12               | 5    |
| MB30235B <sup>7</sup> | 11 x 11 x 5   | 279 x 279 x 127 | Large Bin               | 9                | 4    |
| HK26C                 | 3 <sup>1</sup> / <sub>4</sub> x 2 x 3 <sup>1</sup> / <sub>2</sub> | 95 x 51 x 89    | Double Hook             | .2               | .1   |

<sup>1</sup> Weight capacity = 100 lbs. (45kg) evenly distributed.

<sup>2</sup> 6 outlets, 115V, 15A, 60Hz, illuminated power switch, 15 amp fuse.


<sup>3</sup> UL Listed, accepts 2 - 48" (1219mm), 40W fluorescent bulbs, not included.

<sup>4</sup> Weight capacity = 250 lbs. (113kg) total (including all accessories and bins).

<sup>5</sup> Weight capacity = 50 lbs. (23kg) evenly distributed.


<sup>6</sup> Order in multiples of 12.

<sup>7</sup> Order in multiples of 6.



### Components:

- 1 — Stainless work surface
- 2 — Tie Bar
- 3 — 3-Sided support frame\*
- 4 — Galvanized bottom shelf
- 5 — Foot rest



All Metro Catalog Sheets are available on our website: [www.metro.com](http://www.metro.com)



### InterMetro Industries Corporation

North Washington Street • Wilkes-Barre, PA 18705 • 570.825.2741  
 Fax: 800.638.9263 (East Coast/Canada) • Fax: 800.638.3292 (West Coast)

FOR PRODUCT INFORMATION/CUSTOMER SERVICE:

U.S./Canada/Latin America: 1.800.992.1776 • Europe: +31.76.587.7550

Asia/Pacific: +65.6567.8003 • Middle East/Africa: +971.4.811.8286

Information and specifications are subject to change without notice. Please confirm at time of order

L03-229  
 Rev. 7/12  
 Printed in U.S.A.

Copyright © 2012 InterMetro Industries Corp.

