



Item # _____

Job _____

BASIX PLUS™ — PEDIATRIC CODE RESPONSE CART

All the basics, plus lighter, cleaner and easier.

Basix Plus™ Pediatric Code Response Cart works in conjunction with the Broselow® Pediatric Emergency System. Each cart has 30 inches (762mm) of drawer space, full-extension lightweight aluminum drawers and comes standard with a polymer work surface complete with molded ergonomic handles. Select from a wide range of accessories.

- **Lighter** — Basix Plus carts are made lighter than typical all-steel carts by utilizing aluminum drawers and polymer components such as the polymer base, cart top and corner extrusions.
- **Cleaner** — Polymer corner extrusions provide accessory attachment channels that are accessible for cleaning. The cart's steel shell is not compromised by pre-drilled holes.
- **Easier** — Lighter weight carts are easier to push, and reduce back strain. Accessible accessory channels are easier to clean. Accessory attachment panel provides flexibility for easier accessory positioning.
- **Safe** — The base allows the casters to be mounted outboard of the cabinet providing an anti-tip mechanism. The polymer base doubles as the cart's bumper.
- **Ergonomic** — Push handles are molded into to the polymer cart top providing a natural position to either push or pull the cart. Cart top has a contoured front lip to provide a comfortable work space.
- **Flexible** — User defined configurations can be accommodated by filling an empty cart shell with desired drawer heights. Choose from 3", 6", 9" and 12" (76, 152, 229, 305mm) drawer heights.
- **Accessible** — Full-extension drawers and a full array of accessories put equipment and supplies where needed for convenient access.



MBP8100TN, MB-PEDS
shown with optional accessories



InterMetro Industries Corporation
North Washington Street
Wilkes-Barre, PA 18705
www.metro.com



Basix Plus™ Pediatric Code Response Cart

36.50



BASIX PLUS™ — PEDIATRIC CODE RESPONSE CART

Specifications:

Cart Body:	Double wall construction Inner wall: 22 gauge (.030"/.8mm) steel Outer wall: 20 gauge (.036"/.9mm) steel Finish: Powder coated
Base:	Blow-molded, high-density polyethylene Steel reinforced caster channels Caster mounts located outside of cart body footprint to provide counterbalance to full-extension drawers
Corner Extrusions:	3 channel PVC .125 (3mm) wall thickness
Work Surface:	.125" (3mm) ABS, 19 ¹ / ₂ " x 27 ¹ / ₄ " (495 x 692mm) Molded in handgrips Raised on three sides 1 ⁵ / ₈ " (41mm) Front lip 1/2" (13mm) H Usable work space: 16" x 23" (406 x 584mm) Molded in pencil holder, pill/paper clip holders
Casters:	5" (127mm) diameter with polyurethane tread 2 swivel-brake, 2 swivel-lock (converts from swivel to rigid) 200 lbs. (90kg) per caster, 600 lbs. (270kg) per cart
Drawers:	50 lbs. (23kg) load rating per drawer Body: .032" (.8mm) coated aluminum Finish: Powder coated Fronts: .032" (.8mm) coated aluminum Finish: Powder coated Pull: Polymer Slides: Zinc-Plated Steel Load Rating: 100 lbs. (45kg)/set

All Metro Catalog Sheets are available on our Web Site: www.metro.com



InterMetro Industries Corporation

North Washington Street, Wilkes-Barre, PA 18705
Phone: 570-825-2741 • Fax: 570-825-2852

For Product Information Call: 1-800-433-2232

L03-273
Rev. 02/09

Information and specifications are subject to change
without notice. Please confirm at time of order.

Copyright © 2008 InterMetro Industries Corp.

