



Item # _____

Job _____

BASIX™ — PEDIATRIC CODE RESPONSE CART

Reliable performance for any application.

Basix™ Pediatric Code Response Cart works in conjunction with the Broselow® Pediatric Emergency System. Each cart has 30 inches (762mm) of drawer space, full-extension drawers and comes standard with a polymer work surface and polymer push handles. Select from a wide range of accessories to customize the cart to your needs.

- **Safe** — The polymer base protects the cart and surrounding environment from incidental contact. The base is designed to allow the 5" (127mm) casters to be mounted outboard of the cabinet, which keeps the cart stable when a fully loaded drawer is accessed.
- **Ergonomic** — Polymer top provides a user-friendly work surface. Optional accessory push handles can be mounted in an ergonomic fashion to allow for easy pushing or pulling of the cart.
- **Flexible** — User defined configurations can be accommodated by filling an empty cart shell with desired 3", 6", 9" and 12" (76, 152, 229, 305mm) drawer heights.
- **Accessible** — Full-extension drawers and a full array of accessories put equipment and supplies where needed for convenient access. Optional drawer dividers keep supplies neat and organized.
- **Mobile** — This cabinet comes with 5" (127mm) diameter casters (2 with swivel brakes, 2 with swivel-lock) allowing it to be rolled easily from location to location or moved for convenience when cleaning.
- **Accessories** — A full array of task support accessories are available for above the work surface or the side of the cart.



MBX8100TB, MB-PEDS
shown with optional accessories



InterMetro Industries Corporation
North Washington Street
Wilkes-Barre, PA 18705
www.metro.com



Basix™ Pediatric Code Response Cart

36.60



BASIX™ — PEDIATRIC CODE RESPONSE CART

Specifications:

Cart Body:	Double wall construction Inner wall: 22 gauge (.030"/.8mm) steel Outer wall: 20 gauge (.036"/.9mm) steel Finish: Powder coated
Base:	Blow-molded, high-density polyethylene Steel reinforced caster channels Caster mounts located outside of cart body footprint to provide counterbalance to full-extension drawers
Work Surface:	.125" (3mm) ABS , raised on four sides Dimensions: 11¼" H x 18½" W x 26¾" L (32 H x 470 W x 679mm L) Usable work space: 17" x 25" (432 x 635mm)
Casters:	5" (127mm) diameter with polyurethane tread 2 swivel-brake, 2 swivel-lock (converts from swivel to rigid) 200 lbs. (90kg) per caster, 600 lbs. (270kg) per cart
Drawers:	50 lbs. (23kg) load rating per drawer Body: 22 Gauge (.030"/.8mm) steel with 16 gauge (.060"/2mm) lock clips Finish: Powder coated Fronts: .032" (.8mm) coated aluminum Finish: Powder coated Pull: Polymer Slides: Zinc-Plated Steel Load Rating: 100 lbs. (45kg)/set

All Metro Catalog Sheets are available on our Web Site: www.metro.com



InterMetro Industries Corporation

North Washington Street, Wilkes-Barre, PA 18705
Phone: 570-825-2741 • Fax: 570-825-2852

For Product Information Call: 1-800-433-2232

L03-272
Rev. 02/09

Information and specifications are subject to change
without notice. Please confirm at time of order.

Copyright © 2008 InterMetro Industries Corp.

