




We put space to work.

Item # _____

Job _____

Lionville® Series Medication Carts

- The lightweight design of the carts make them ideal for medication transport and distribution.
- The range in sizes enables the cart system to be configured for a few patients for bedside application to high capacity to serve as the main medication storage hub.
- The flexible design allows for a number of user defined options that can be integrated into the cart including med cup dispenser storage wells and bins, additional storage, or concealed trash.
- Security is key for medication storage and the Lionville series offers a number of options starting with a basic key lock that can be upgraded to include:
 - Autolock supporting basic key function
 - Electronic keyless entry
 - Electronic keyless entry with autolock
 - Electronic keyless entry, LCD screen, autolock and audit capabilities*
 - Keyed double lock for narcotics
 - Keyless double lock for narcotics*
- Carts can be configured with a mixture of cassettes, storage drawers or narcotics drawers to suit specific needs.
- Drawer pulls are available in five colors to complement décor or differentiate the carts.



All keyless entry options include keyed access as a backup for entry.

*Keyless narcotic storage compatible with LCD screen keyless entry version only.



800



600



400

All Metro Catalog Sheets are available on our website: www.metro.com



InterMetro Industries Corporation
 North Washington Street, Wilkes-Barre, PA 18705
 Product Information. U.S. and Canada: 1.800.992.1776
 Outside U.S. and Canada: www.metro.com/contactus

L02-192e
Printed in U.S.A. Rev. 9/17

Information and specifications are subject to change without notice. Please confirm at time of order.

Copyright © 2017 InterMetro Industries Corp.



Lionville® Series Medication Carts

37.30

Lionville® Series Medication Carts




We put space to work.®

Features:

- Lightweight aluminum construction
- Scratch resistant powder-coat finish
- Polymer molded top
- Wrap-around base bumper and full height corner bumpers
- 5" (127mm) diameter casters, 2 total lock (F), 2 swivel lock (R)
- Front control panel with mechanical main lock featuring the patented* drawer grabber latching mechanism, key removable core and color-coded re-lock handle
- Optional double-lock drawers, keyless entry, automatic re-lock and reporting
- Integrated side storage recess or optional pullout work surface
- Integrated rear compartment for storage, trash or technology integration



All carts are 23½" (597mm) deep x 41⅞" (1045mm) high.



	LMED800/LMED800R	LMED600/LMED600-27	LMED400/LMED400R
Width	22" (559mm)	32⅝" (829mm)	39⅞" (994mm)
Capacity	9 tiers	10 tiers** 28" (711mm)	19 tiers
Casters, options	5" (127mm) diameter	5" (127mm) diameter	5" (127mm) diameter
Storage Drawers 11⅝" x 15⅞" (295x384mm)	3" (76mm) H, 6" (152mm) H, 9" (229mm) H, 12" (305mm) H		
Side Drawers 600 19¼" x 8⅝" (489x219mm)		4" (102mm), 8" (203mm), 10" (254mm)	
Cassettes	13"D x 16" W (330 x 406mm) One Tier: 3⅞" H (79mm) Two Tier: 6⅝" H (161mm)		
Patient Drawers	11⅞" D x 2½" H 3 per tier: 5" W (127mm), 2 per tier: 7½" W (191mm)		

*US Patent # 6511138

**Optional 18" w (457mm) pullout work surface will default a 600 cart to 9 tier.

Options:

- Power and computer integration
- Pullout side shelves (600 & 400 only)
- Pullout side work surface
- Narcotics drawer
- Integrated or side waste bins
- Worksurface well or drop in utility bin
- Sharp container bracket
- Tilt-out side bins
- Dividers
- Light
- Cup dispensers
- Glove box holder

an Ali Group Company



The Spirit of Excellence